

Chemblend² 2-Part Developer

***Balanced Speed
and Contrast
with Low Fog
for Maximum
Diagnostic Detail
on All X-ray Film***



If you have enjoyed the detail that Chemblend[®] 3-part Developer gave you, but didn't enjoy the shipping surcharge for the acetic acid, now you can have it both ways!

H.R. Simon has kept the working formula the same, but now **Chemblend² 2-part Developer** is DOT non-hazardous, which eliminates UPS HAZMAT surcharges.

Chemblend² 2-part Developer has an excellent shelf life.

Radiology departments consistently report excellent Quality Assurance, equalling the top performance of Chemblend 3-part developer.



It's fair to say that **Chemblend² 2-part Developer** is the new industry standard for liquid concentrate x-ray developer. You won't have to change a thing about how you process – you'll just have to get used to not paying more for shipping.

CBM20
Chemical
Mixing
Station



www.hrsimon.com

H.R. Simon and Company, Inc.
3515 Marmenco Court
Baltimore, Maryland 21230 • USA

TOLL FREE 800-638-9460
BALTIMORE 410-636-5555
FAX 410-636-4415
E-MAIL sales@hrsimon.com

Chemblend² 10 Gallon Carton
Stock Number

S60-045

*2 bottles, 5 quart part A developer
2 bottles, 16 ounce part B developer*

S60-045		CB ² Developer 10 Gallon Pack			31 lbs / 14.1 kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS	
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B
5	19	3	80	13.7	1	1
10	38	7	32	27.5	2	2
S50-094		CB ² Developer 200 Gallon Bulk Mix			609 lbs / 277 kg	
TO MAKE		USE THIS AMOUNT OF WATER			USE THESE UNITS	
Gallons	Liters	Gallons +	Ounces	Liters	Part A	Part B
200	757	145	0	548.8	1	1

*Never leave containers half-full

PRODUCTS • Chemblend Developers • SEF² Developer • Autochem[®] Dental • HS Manual Solutions
• CBM20 Mixing Station • Imaging Services • Equipment and Accessories